

overbed tables

Amico Furniture is pleased to announce the addition of Overbed tables to its line of healthcare products. With our unique quality and design features, Amico's Overbed tables allow the end user to fully customize the product to fulfill every need.



At Amico Furniture, it is our goal to provide a variety of product features to complement every décor.

Features

- 18" x 32" Table top slides forwards or backwards to reveal lower 17" x 26" surface
- Easily accessible tabs for table top slide activation by patient or caregiver
- Adjust height easily with automatic gas cylinder lift mechanism
- Non breakable acrylic mirror
- Laminate top resist damages from spills and overall use
- 2" twin wheel casters make for smooth easy handling
- Supports 100 lbs. distributed evenly on top table surface

Options

- Choice of frame color
- Customize with 9 different laminate selections
- U Base or H base option
- Spill guard to ensure containment of liquid
- Lower base (clearance 2.5")
- Protective guard around edge of mirror available



Base Selections



Laminate Selections



the right connection