

## CONVERSION TABLE FOR WALL MOUNTED OUTLETS

	OHIO MEDICAL/OHMEDA		PURITAN BENNETT			NCG/CHEMETRON		OXEQIP/MEDSTAR	
	DIAMOND 1	DIAMOND 2 & 3	SERIES 73	SERIES "A"	LITHARGE	SERIES 148,400,460	SERIES 500, 560	SERIES # 7 OR 5/8-24	SMOOTH BORE OR 5/8-27
INSTALLATION TOOL	O-VC-OH1-TOOL	Not needed	O-CPBTOOLKIT	Not needed	Not needed	O-VC-OH1-TOOL	Not needed	O-VC-OH1-TOOL	Not needed
KIT NUMBERS		<b>Complete Conversion</b>							
QD	O-OH1-QD-L-G	O-OH2-QD-L-G	O-C73-QD-L-G-W	O-CSA-QD-L-G-W	O-CLI-QD-L-G-W	O-NCG-QD-L-G-W	O-NCG500-QD-L-G	O-OX1-QD-L-G	O-SBCO-QD-L-G
DI	O-OH1-DI-L-G	O-OH2-DI-L-G	O-C73-DI-L-G-W	O-CSA-DI-L-G-W	O-CLI-DI-L-G-W	O-NCG-DI-L-G-W	O-NCG500-DI-L-G	O-OX1-DI-L-G	O-SBCO-DI-L-G
CH	Not Available	Not Available	Not Available	Not Available	Not Available	Not Available	Not Available	O-OX1-CH-L-G	Not Available
PB	Not Available	Not Available	O-C73-PB-L-G-W	O-CSA-PB-L-G-W	O-CLI-PB-L-G-W	O-NCG-PB-L-G-W	O-NCG500-PB-L-G	O-OX1-PB-L-G	O-SBCO-PB-L-G
		<b>Using Existing Parts</b>							
		O-OH2-QD-L-G-PP							
		O-OH2-DI-L-G-PP							

**L** defines the Language:

English (NFPA) = U  
 English (CSA) = E  
 French (CSA) = F  
 Spanish (NFPA) = S

**G** defines the Gas Type:

Oxygen = O  
 Medical Air = A  
 Nitrous Oxide = 2  
 Carbon Dioxide\* = C  
 Nitrogen\*\* = N  
 MedVac = V  
 Waste Anesthetic Gas Disposal = W (NFPA)  
 Anesthetic Gas Scavenging System = AG (CSA)

**W** defines the Frame Width:

3 = 3.25"  
 5 = 5.0"

\* Only Available in DISS or QD

\*\* Only Available in DISS

**VERSION 1:12**