



reflection series

artwall


product benefits

- conceals medical gas outlets and equipment for a non-clinical appearance
- over 500 prints and frames offered to coordinate and enhance room décor
- swiveling outlets allow for frame to close while in use
- diagonal or horizontal movement to accommodate your room design

applications


LDR-Post-Partum
Patient Rooms
Short Stay
Hospice
Oncology

product features



- ① a wide selection of art prints and frames available
- ② smooth movement in horizontal or diagonal, left hand or right hand direction
- ③ customizable back box depth
- ④ Amico's swiveling outlet allows for easy access during standard maintenance
- ⑤ available 2,3,4 or 5 gas configurations
- ⑥ rough-in assembly is easily installed in the stud work of the wall
- ⑦ the diagonal movement ArtWall includes a safety piston to prevent accidental closure on fingers

typical drawing



Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more.

represented by: