



# regal series

## surface mount flatwalls


### product benefits

- vertical configuration to reduce equipment overhang
- increased equipment management with four integrated vertical accessory rails per unit
- greater flexibility allowing services on all sides of the unit
- ease of installation with medical gas and electrical services terminated to a single point at or above the ceiling line
- expandable to meet your future requirements

### applications


ICU  
NICU  
Emergency Department  
LDR Post-Partum  
PACU  
Transitional Step-Down  
Patient Rooms  
Radiology  
Special Procedure

**product features**



- ① Available in full, 3/4, 5/8, and 1/2 length options
- ② 19", 28" and 38" widths available. Custom widths can be manufactured to meet your needs
- ③ four integrated, scratch resistant vertical rails included
- ④ outlet and electrical receptacle layout based on facility needs
- ⑤ overbed light and bed locator can be integrated to the unit through a surface mounted raceway
- ⑥ over 2000 colors of hospital grade laminate to select from including Formica, Wilson Art, Nevamar and Pionite
- ⑦ factory pre-wired and piped to a single point connection for ease of installation

**typical drawing**



Please visit [www.amico.com](http://www.amico.com) for all technical specifications, installation instructions, detailed AutoCAD drawings and more.

**represented by:**