

overbed light



bed locator



overbed light


- ▶ micro-switch bed-stop safety feature to prevent damage from rising bed
- ▶ optional low voltage controller provides easily accessible switching
- ▶ over 2000 colors of hospital grade laminate to select from Formica, Wilson Art, Nevamar and Pionite to meet the rooms décor and match headwall system
- ▶ overbed and ceiling mount options available in 3 and 5 function

bed locator

- ▶ single piece high impact resistant Kydex to increase durability
- ▶ conduit raceways allow for custom configurations and future modifications of electrical and communication services
- ▶ designed to safely position and secure hospital bed preventing damage to walls

overbed lights and bed locators


product features




- ① sound rated ballasts for separate up and down light operation
- ② high impact resistant clear ribbed lens
- ③ pull down handles swivel light downwards toward patient area for bedside exam
- ④ automatic bed-stop interlock safety feature

product features

- ① made from single-piece high impact vacuum formed Kydex material
- ② accommodates up to four electrical/communication devices on each side of the unit
- ③ conduit raceways used for low voltage, normal and emergency power



typical drawing



Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more.

represented by: