



noble series

neonatal intensive care unit (NICU)


product benefits

- fully customizable to meet caregiver's needs
- vertical and horizontal accessory rails for increased equipment management
- adaptable to level 1, 2, and 3 Neonatal Intensive Care Units
- cabinet work can be included for additional counter, shelving, and storage space
- removable panels for maintenance and future expansion
- medical gas and electrical services positioned at various angles and heights for unobstructed accessibility

applications


NICU
(Neonatal Intensive Care Unit)

product features



- ① available in configurations such as back to back, stand alone and against the wall
- ② scratch resistant, integrated accessory rails come in horizontal and vertical configurations
- ③ factory pre-wired and piped to a single-point connection for ease of installation
- ④ cabinet work customized to meet your specific requirements
- ⑤ can accommodate a large quantity of medical gas and electrical services in various configurations
- ⑥ includes a removable 3-drawer cart
- ⑦ over 2000 colors of hospital grade laminate to select from including Formica, Wilson Art, Nevamar and Pionite

typical drawing



Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more.

represented by: