

Amico Furniture is pleased to announce the addition of Recliner and Treatment chairs to its line of healthcare products. With our unique quality and design features, Amico’s Recliner and Treatment chairs allow the end user to fully customize the product to fulfill every need. At Amico Furniture, it is our goal to provide a variety of product features to compliment every décor.



recliner chairs

Features

- 3 position recline, including full layout position
- Safer small opening between foot rest and seat
- Great foot rest stability
- Maple wood arm caps available
- Optional IV rod and flip-up tablet

treatment chairs

Features

- Steel framed chairs with molded contour backrest
- Standard models measures 20" between arms
- Extra wide 24" model available
- seat cushion and footrest cushion is one piece with no spaces between leg and foot rest
- 3" heavy-duty non-marking casters, rear casters are locking